# ADVANCIS

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
## Actilite®

Non-adherent viscose net dressing coated with 99%  
Manuka honey and 1% Manuka oil

### Product description

Actilite® is an effective broad spectrum antimicrobial dressing containing 99% Manuka honey and 1% Manuka oil for use on low to moderate exuding wounds. The dressing is designed to protect a wound, promote healing and allow the passage of exudate. It is an ideal choice for epithelialising wounds that are perceived to be at risk of re-infection.

### Indications

May be used on all wound types including: Superficial wounds, partial or full thickness wounds, e.g. cuts and abrasions, pressure ulcers, leg ulcers, surgical wounds, burns, infected wounds, graft & donor sites. May assist with debridement if used in conjunction with Activon Tube.

### Directions for use

Remove the clear film backing from both sides of the dressing. Place directly onto wound bed either side down, may be cut to size. Cover with a secondary dressing, depending on tissue type within the wound bed and level of exudate.

### Frequency of change

Depending on wound exudate levels, any surrounding interstitial fluid, oedema and dressing regime the dressing may be left in place for up to 7 days. Will be less effective when diluted by wound exudate. To maintain an effective honey level, change the dressing when the colour changes significantly. Facilitate removal of the dressing by irrigating the wound bed with saline solution.

### Warning

Device is intended as a single-use device. Re-using the device presents the risk of infection and is not recommended.

### Contra-indications

Do not use if patient is allergic to bee venom, bee products or essential oils. Actilite® contains glucose. Monitor the levels of patients with diabetes. A few patients experience pain when applied, depending on sensitivity of the wound it may be necessary to consider an appropriate level of analgesia. If pain continues unabated, discontinue use and irrigate the wound with saline solution.

### Warnings and Precautions

Actilite® is a single use product. If used on more than one patient, cross-contamination or infection may result. Opening the dressing pack compromises the sterile barrier. Any unused dressings should be discarded.

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