ADVANCIS

ADVANCIS

Brightwake Ltd, Lowmoor Business Park, Kirkby in Ashfield, Nottinghamshire, NG17 7JZ
Tel: +44 (0) 1623 751500
Fax: +44 (0) 871 264 8238
Email: info@advancis.co.uk
Web: www.advancis.co.uk

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
ADVANCIS

ADVANCIS

	Catalogue Number		Consult instructions for use
	Do not re-use		Keep dry
	Sterilized using irradiation		Keep away from sunlight
	Do not use if package is damaged		Manufacturer
	Caution		Authorized representative in the European Community
	Quantity		Batch code
	Unique Device Identifier (UDI) carrier including the AIDC and human readable information)		Medical Device
	Date of manufacture/ Country of manufacture		CE Mark for Northern Ireland, EC & EU
	The product has made a financial contribution towards the recovery and recycling of packaging in Europe		Importer
	Do not re-sterilize		Single sterile barrier system
	Use-by date		Indicates the authorised representative in Switzerland

Activon Tulle

Knitted viscose mesh impregnated with 100% Manuka honey

Product description

Activon® Tulle is a knitted viscose mesh primary dressing impregnated with 100% Medical Grade Manuka honey. The dressing is designed to protect a wound, promote healing and allow the passage of exudate. It is an ideal choice for epithelialising wounds that are perceived to be at risk of re-infection.

Indications

May be used on all wound types including: Superficial wounds, partial or full thickness wounds, e.g. cuts and abrasions, pressure ulcers, leg ulcers, surgical wounds, burns, infected wounds, graft & donor sites. May assist with debridement if used in conjunction with Activon Tube.

Directions for use

Remove the clear film backing from both sides of the dressing. Place directly onto wound bed either side down, may be cut to size. Cover with a secondary dressing, depending on tissue type within the wound bed and level of exudate.

Frequency of change

Depending on wound exudate levels, any surrounding interstitial fluid, oedema and dressing regime the dressing may be left in place for up to 7 days. Will be less effective when diluted by wound exudate. To maintain an effective honey level, change the dressing when the colour changes significantly. Facilitate removal of the dressing by irrigating the wound bed with saline solution.

Warning

Device is intended as a single-use device. Re-using the device presents the risk of infection and is not recommended.

Contra-indications

Do not use if patient is allergic to bee venom, bee products or essential oils. Activon® Tulle contains glucose. Monitor the levels of patients with diabetes. A few patients experience pain when applied, depending on sensitivity of the wound it may be necessary to consider an appropriate level of analgesia. If pain continues unabated, discontinue use and irrigate the wound with saline solution.

Warnings and Precautions

Activon® Tulle is a single use product. If used on more than one patient, cross-contamination or infection may result. Opening the dressing pack compromises the sterile barrier. Any unused dressings should be discarded.

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