Scan for more information

**Adv**ancis  
Medical

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Medical

## Advazorb®/Advazorb® Lite

Hydrophilic foam dressing with film backing

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
pt BR

### Product description

A soft, conforming non-adherent hydrophilic polyurethane foam wound dressing with pink polyurethane film backing. Advazorb® has high fluid absorbency and retention characteristics.

### Indications

Suitable for a wide range of exuding wounds, granulating wounds, leg ulcers, pressure sores and for protection of bony prominences.

### Contra-indications

Arterial bleeds and heavily bleeding wounds.

### Directions for use

Place directly to the wound surface pink film side up, ensuring the absorbent pad covers the entire wound area with an overlap of 1-2cm. Secure in place with tape, appropriate bandage or film dressing. Can be used under compression.

### Frequency of change

Can be left in place for up to 7 days but should be changed immediately if strike through is apparent or when dressing reaches capacity. If exudate is visible around the edges of the dressing this is a clear indication that a dressing change is required. Clinical observation is necessary to determine required frequency of change.

### Warnings and Precautions

Advazorb® is a single use product. If used on more than one patient, cross-contamination or infection may result. Opening the dressing pack compromises the sterile barrier. Any unused dressings should be discarded.

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