ADVANCIS

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Superabsorbent foot wrap dressing

Product description

Eclypse® Foot is an exudate management wound dressing especially designed for the foot. The dressing combines a rapid wicking layer, a CrystaLock™ layer, which is highly absorbent and locks away exudate, a sub-backing layer for even exudate distribution and a non-strikethrough backing to prevent leakage.

Indications

Moderate to heavily exuding wounds such as leg ulcers, pressure ulcers and diabetic foot ulcers.

Contra-indications

Arterial bleeds and heavy bleeding wounds.

Directions for use

Open out the dressing fully, white side up placing the end of foot in the middle of the dressing so the heel is at the edge. Fold the side panels over the foot one at a time and secure with tape. Do not cut Eclypse® Foot. Eclypse® Foot can be used under compression bandaging. Eclypse® Foot can be used as a primary or secondary dressing.

Frequency of change

Wear time will depend on the level of exudate. Eclypse® Foot must be changed before dressing reaches capacity but depending on wound exudate levels, may be left in place for up to 7 days.

Disposal

Dispose as clinical waste.

Do not re-sterilize

Single sterile barrier system

Use-by date

Indicates the authorised representative in Switzerland

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