MEDICA PACKAGING

CAD SPECIFICATION

GENERAL INFORMATION

CAD SPEC NO	P028522
CUSTOMER	Brightwake
ACCOUNT NO	BRIGHTWAKE
CUSTOMER REF	

DESIGN INFORMATION

MATERIAL	Incada Silk 540mic
CALIPER	0.54mm
GRAMMAGE (gsm)	325
STYLE	RTE
LENGTH	295mm
WIDTH	35mm
HEIGHT	330mm
BLANK WIDTH	672mm
BLANK HEIGHT	438mm
BLANK AREA	256,181.66sq mm
BLANK WEIGHT	**83.26gm**

CAD DESIGN NOTES

VERSION INFORMATION

CREATED BY	DAVE
CREATED	01/05/2026
MODIFIED BY	DAVE
MODIFIED	01/05/2026
VERSION	1

REASON FOR LATEST VERSION

CAD LINE TYPE LEGEND

CUT		BRAILLE	
CREASE		ACCUBRAILLE	
SCORE		EMBOSS	
REV SCORE		VARNISH FREE	
PERF		TECHNICAL	

Approved by C.Allsop on 19th May

ADVANCIS
ADVAZORB BORDER

Hydrophilic foam dressing
with soft silicone contact layer and border
20cmx20cm (8inx8in) **x10**

ADVANCIS
ADVAZORB BORDER

Hydrophilic foam dressing
with soft silicone contact layer and border
20cmx20cm (8inx8in) **x10**

X10

ADVANCIS
ADVAZORB BORDER

Hydrophilic foam dressing
with soft silicone contact layer
and border
20cmx20cm (8inx8in)

REF **CR4195**

Silfix TECHNOLOGY

AW199/v9

ADVBOX4195/10

Hydrophilic foam dressing
with soft silicone contact layer and border
20cmx20cm (8inx8in) **x10**

ADVANCIS
ADVAZORB BORDER

Hydrophilic foam dressing
with soft silicone contact layer and border
20cmx20cm (8inx8in) **x10**

ADVANCIS
ADVAZORB BORDER