MEDICA PACKAGING

CAD SPECIFICATION

GENERAL INFORMATION

CAD SPEC NO	P028503
CUSTOMER	Brightwake
ACCOUNT NO	BRIGHTWAKE
CUSTOMER REF	

DESIGN INFORMATION

MATERIAL	Incada Silk 540mic
CALIPER	0.54mm
GRAMMAGE (gsm)	325
STYLE	RTE
LENGTH	270mm
WIDTH	50mm
HEIGHT	280mm
BLANK WIDTH	652mm
BLANK HEIGHT	428mm
BLANK AREA	227,976.89sq mm
BLANK WEIGHT	74.09gm

CAD DESIGN NOTES

VERSION INFORMATION

CREATED BY	DAVE
CREATED	01/05/2026
MODIFIED BY	DAVE
MODIFIED	01/05/2026
VERSION	1

REASON FOR LATEST VERSION

CAD LINE TYPE LEGEND

CUT		BRAILLE	
CREASE		ACCUBRAILLE	
SCORE		EMBOSS	
REV SCORE		VARNISH FREE	
PERF		TECHNICAL	

Claire Allsop approved 26th May 2026

ADVANCIS
ECLYPSE ADHERENT

Super Absorbent Dressing with Silicone Contact Layer
15cmx15cm (6inx6in) **X10**

ADVANCIS ECLYPSE ADHERENT

Super Absorbent Dressing with Silicone Contact Layer
15cmx15cm (6inx6in) **X10**

X10

ADVANCIS
ECLYPSE ADHERENT

Super Absorbent Dressing with Soft Silicone Contact Layer
15cmx15cm (6inx6in)

REF **CR3863**

Silfix® TECHNOLOGY

Super Absorbent Dressing with Silicone Contact Layer
15cmx15cm (6inx6in) **X10**

ADVANCIS ECLYPSE ADHERENT

AW016/v10

ADVBOX3863/11