MEDICA PACKAGING

CAD SPECIFICATION

GENERAL INFORMATION

CAD SPEC NO	P028515
CUSTOMER	Brightwake
ACCOUNT NO	BRIGHTWAKE
CUSTOMER REF	

DESIGN INFORMATION

MATERIAL	Incada Silk 540mic
CALIPER	0.54mm
GRAMMAGE (gsm)	325
STYLE	RTE
LENGTH	296mm
WIDTH	30mm
HEIGHT	350mm
BLANK WIDTH	664mm
BLANK HEIGHT	448mm
BLANK AREA	263,249sq mm
BLANK WEIGHT	85.56gm

CAD DESIGN NOTES

VERSION INFORMATION

CREATED BY	DAVE
CREATED	01/05/2026
MODIFIED BY	DAVE
MODIFIED	01/05/2026
VERSION	1

REASON FOR LATEST VERSION

CAD LINE TYPE LEGEND

CUT		BRAILLE	
CREASE		ACCUBRAILLE	
SCORE		EMBOSS	
REV SCORE		VARNISH FREE	
PERF		TECHNICAL	

Approved by C.Allsop on 19th May

ADVANCIS
ECLYPSE ADHERENT SACRAL

Superabsorbent sacral dressing
with Silfix soft silicone contact layer
22cmx23cm (8.5inx9in) **x10**

ADVANCIS
ECLYPSE ADHERENT SACRAL

Superabsorbent sacral dressing
with Silfix soft silicone contact layer
22cmx23cm (8.5inx9in) **x10**

X10

ADVANCIS
ECLYPSE ADHERENT SACRAL

Superabsorbent sacral dressing
with Silfix soft silicone contact layer
22cmx23cm (8.5inx9in)



Superabsorbent sacral dressing
with Silfix soft silicone contact layer
22cmx23cm (8.5inx9in) **x10**

REF **CR3986**

CrystaLock TECHNOLOGY

ADVANCIS
ECLYPSE ADHERENT SACRAL

AW024/v10

ADVBOX3986/10